

Product Datasheet

Human SARM1 protein (orb595102)

Description

Recombinant Human Sterile alpha and TIR motif-containing protein

Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Full Length of Mature Protein
Protein Sequence	LAVPGPDGGGGTGPWWAAGGRGPREVSPGAGTEVQDALERALPELQQALSAL KQAGGARAVGAGLAEVFQLVEEAWLLPAVGREVAQGLCDAIRLDGGLDLLRL QAPLETRVQAARLLEQILVAENRDRVARIGLVILNLAKEREPVELARSVAGILE HMFKHSEETCQRLVAAGGLDAVLYWCRRTDPALLRHICALALGNALHGGQAV QRRMVEKRAAEWLFPLAFSKEDELLRLHACLAVAVLATNKEVEREVERSGTLAL VEPLVASLDPGRFARCLVDASDTSQGRGPDDLQRLVPLLDNRLEAQCIGAFYL CAEAAIKSLQGKTKVFSDIGAIQSLKRLVSYSTNGTKSALAKRALRLLGEEVPRPIL PSVPSWKEAEVQTLWQQIGFSKYCESFREQQVDGDLRLTEELQTDLGMKSG ITRKRFFRELTELKTFANYSTCDRNLADWLGLDPRFRQYTYGLVSCGLDRSLL HRVSEQQLLEDGCIHLGVHRARILTAAREMLHSPLCTGGKPSGDTDPVFISYRR NSGSQLASLLKVHLQLHGFSVFIDVEKLEAGKFEDKLIQSVMGARNFVLVSPGA LDKCMQDHDCKDWVHKEIVTALSCGKNIVPIIDGFEWPEPQVLPEDMQAVLTFN GIKWSHEYQEATIEKIRFLQGRSSRDSSAGSDTSLEGAAPMGPT
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Baculovirus
MW	80.4 kDa
Uniprot ID	Q6SZW1
Expiration Date	6 months from date of receipt.