

Product Datasheet

Human GLP-1 protein (orb595165)

Description	Recombinant Human Glucagon-like peptide 1 receptor(GLP1R)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged
Note	For research use only
Application notes	Full Length of Mature Protein
Protein Sequence	RPQGATVSLWETVQKWREYRRQCQRSLTEDPPPATDLFCNRTFDEYACWPDG EPGSFVNVSCPWYLPWASSVPQGHVYRFCTAEGWLQKDNSSLPWRLSECEE SKRGERSSPEEQLLFLYIIYTVGYALSFSALVIASAILLGFRHLHCTRNYIHLNLFASF ILRALSVFIKDAALKWMYSTAAQQHQWDGLLSYQDSLSCRLVFLLMQYCVAANY YWLLVEGVLYTLAFSVLSEQWIFRLYVSIGWGVPLLFVVPWGIVKLYEDEGC WTRNSNMNYWLIIRLPILFAIGVNFLIFVRVICIVVSKLKANLMCKTDIKCRLAKSTL TLIPLLGTHEVIFAFVMDEHARGTLRFIKLFTLSFTSFQGLMVAILYCFVNNEVQL EFRKSWERWRLEHLHIQRDSSMKPLKCPTSSLSSGATAGSSMYTATCQASCS
Purity	Greater than 85% as determined by SDS-PAGE.
Source	in vitro E.coli expression system
MW	53.6 kDa
Uniprot ID	P43220
Expiration Date	6 months from date of receipt.