

Product Datasheet

E. coli ompA protein (orb595197)

Description	Recombinant Escherichia coli Outer membrane protein A(ompA)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged
Note	For research use only
Application notes	Full Length of Mature Protein
Protein Sequence	APKDNTWYTGAKLGWSQYHDTGFINNNGPTHEQLGAGAFGGYQVNPYVGFE MGYDWLGRMPYKGSVENGAYKAQGVQLTAKLGYPITDDLDIYTRLGGMVWRA DTKSNVYGKNHDTGVSPVFAGGVEYAITPEIATRLEYQWTNNIGDAHTIGTRPDN GMLSLGVSYRFGQGEAAPVVAPAPAPAPEVQTKHFTLKSDVLFNFKATLKPEG QAALDQLYSQLSNLDPKDGSVVVLGYTDRIGSDAYNQGLSERRAQSVVDYLISK GIPADKISARGMGESNPVTGNTCDNVKQRAALIDCLAPDRRVEIEVKGIKDVVTQ PQA
Purity	Greater than 85% as determined by SDS-PAGE.
Source	in vitro E.coli expression system
MW	38.0 kDa
Uniprot ID	P0A910
Expiration Date	6 months from date of receipt.