

Product Datasheet

DARPP32 (phospho-Thr75) antibody (orb5961)

Description	DARPP32 (phospho-Thr75) antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human DARPP32 around the phosphorylation site of Thr75 (AY(p-T)PP)
Target	Phospho-DARPP32 (Thr75)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32 kDa
Uniprot ID	Q9UD71
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.