

Product Datasheet

UBE2F Immunizing peptide (orb59728)

Description	UBE2F Immunizing peptide
Preservatives	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Storage	Shipped at ambient temperature, store at - 20°C
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb59728 (NEDD8 conjugating enzyme (N terminus) antibody)
Protein Sequence	SKLKRDDGLKGSRT
Expiration Date	12 months from date of receipt.