


Product Datasheet

eIF2 alpha (Phospho-Ser51) antibody (orb5998)

Description	Rabbit polyclonal antibody to eIF2 alpha (Phospho-Ser51) which belongs to the eIF-2-alpha family. El...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 whole cell...