

Product Datasheet

eIF4G (phospho-Ser1108) antibody (orb6000)

Description	eIF4G (phospho-Ser1108) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human eIF4G around the phosphorylation site of ser1185 (KR(p-S)FS)
Target	phospho-eIF4G (ser1185)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	176 kDa
Uniprot ID	Q04637
Expiration Date	12 months from date of receipt.