

Product Datasheet

Human S100A11 protein (orb6009)

Description	Recombinant of human S100A11 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, 1mM DTT, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Protein Sequence	Met1-Thr105
Source	E.coli
Uniprot ID	P31949
Expiration Date	6 months from date of receipt.