

Product Datasheet

ETK (phospho-Tyr40) antibody (orb6025)

Description	ETK (phospho-Tyr40) antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human ETK around the phosphorylation site of Tyr40 (LS(p-Y)YE)
Target	Phospho-BMX (Tyr40)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	78 kDa
Uniprot ID	P51813
Expiration Date	12 months from date of receipt.