

## Product Datasheet

### POLDIP3 (orb602634)

|                           |                                 |
|---------------------------|---------------------------------|
| <b>Description</b>        | POLDIP3                         |
| <b>Form/Appearance</b>    | 10 µg plasmid + 200µl Glycerol  |
| <b>Note</b>               | For research use only           |
| <b>Uniprot ID</b>         | BC019643                        |
| <b>Hazard Information</b> | Non-Toxic                       |
| <b>Expiration Date</b>    | 12 months from date of receipt. |