

Product Datasheet

FADD (phospho-Ser191) Mouse Specific antibody (orb6030)

Description	FADD (phospho-Ser191) Mouse Specific antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from mouse FADD around the phosphorylation site of Ser191 (NM(p-S)PV)
Target	Phospho-FADD (Ser191)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	23 kDa
Uniprot ID	Q13158
Expiration Date	12 months from date of receipt.