

Product Datasheet

Human C8G protein (orb603949)

Description	Recombinant Human Complement component C8 gamma chain(C8G)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Full length
Protein Sequence	QKPQRPRRPASPISTIQPKANFDAQQFAGTWLLVAVGSACRFLQEQQGHRAEATT LHVAPQGTAMAVSTFRKLDGICWQVRQLYGDTGVLGRFLLQARDARGAVHVVV AETDYQSFAVLYLERAGQLSVKLYARSLPVSDSVLSGFQQRVQEAHLTEDQIFYF PKYGFCEAADQFHVLDDEVRR
Purity	Greater than 85% as determined by SDS-PAGE.
Source	E.coli
MW	26.4 kDa
Uniprot ID	P07360
Expiration Date	6 months from date of receipt.