

Product Datasheet

Virus CYS4 protein (orb603959)

Description	Recombinant <i>Saccharomyces cerevisiae</i> Cystathionine beta-synthase protein(CYS4)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Full Length
Protein Sequence	MTKSEQQADSRHNVIDLVGNTPLIALKKLPKALGIKPQIYAKLELYNPGGSIKDRIA KSMVEEAEASGRIHPSRSTLIEPTSGNTGIGLALIGAIGYRTIITLPEKMSNEKVS LKGALGAEIIRTPAAAWDSPESHIGVAKKLEKEIPGAVILDQYNNMMNPEAHYFGT GREIQRQLEDNLFDNLRAVVAGAGTGGTISGISKYLKEQNDKIQIVGADPFGSIL AQPENLNKTDITDYKVEGIGYDFVPQVLDRLIDVWYKTDDKPSFKYARQLISNE GVLVGSSGSAFTAVVKYCEDHPELTEDDVIVAIFPDSIRSYLTKFVDDDEWLKKN NLWDDDDVLARFDSSKLEASTTKYADVFGNATVKDLHLKPVVSVKETAKVTDVIKI LKDNGFDQLPVLTEGKLSGLVTLSELLRKLSINNSNNDNTIKGKYLDFKKLN NDVSSYNENKSGKKKFIKFDENSKLSDLNRRFEKNSSAVITDGLKPIHIVTKMDLL SYLA
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	60 kDa
Uniprot ID	P32582
Expiration Date	6 months from date of receipt.