

Product Datasheet

Human CYP2C8 protein (orb604035)

Description	Recombinant Human Cytochrome P450 2C8(CYP2C8)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Full Length of Mature Protein
Protein Sequence	RKLPPGPTPLPIIGNMLQIDVKDICKSFTNFSKVYGPVFTVYFGMNPVVFHGYEAV KEALIDNGEEFSGRGNSPISQRITKGLGISSNGKRWKEIRRFSLTLRNFGMGKR SIEDRVQEEAHCLVEELRKTAKSPCDPTFILGCAPCNVICSVVVFQKRFDYKDQNFL TLMKRFNENFRILNSPWIQVCNNFLLIDCFPGTHNKVLKNVALTRSYIREKVKEH QASLDVNNPRDFIDCFLIKMEQEKDNQKSEFNENLVGTVADLFVAGTETTSTTL RYGLLLLLKHPEVTAKVQEEIDHVIGRHRSPCMQDRSHMPYTDVAVVHEIQRYSDL VPTGVPHAVTTDTKFRNYLIPKGTTIMALLTSVLHDDKEFPNPNIFDPGHFLDKNG NFKKSDYFMPFSAGKRICAGEGLARMELFLTTILQNFNLKSVDDLKLNLTAVT KGIVSLPPSYQICFIPV
Purity	Greater than 85% as determined by SDS-PAGE.
Source	E.coli
MW	60.1 kDa
Uniprot ID	P10632
Expiration Date	6 months from date of receipt.