

Product Datasheet

Human ERG protein (orb604077)

Description	Recombinant Human Transcriptional regulator ERG(ERG)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Full Length of Isoform 5
Protein Sequence	MASTIKEALSVMSEDQSLFECAYGTPHLAKTEMTASSSSDYGQTSKMSRPVPQQ DWLSQPPARVTIKMECNPSQVNGSRNSPDECSVAKGGKMGSPDTVGMNYGS YMEEKHMPPPNMTTNERRVIVPADPTLWSTDHVRQWLEWAVKEYGLPDVNILL FQNI DGKELCKMTKDDFQRLTPSYNADILLSHLHYLRETPLPHLTSDDVDKALQN SPRLMHARNTGGAAFI PNTSVYPEATQRITTRPDLPYEP PRRSAWTGHGHPTPQ SKAAQPSPSTVPKTEDQRPQLDPYQILGPTSSRLANPGSGQIQLWQFLLELLSDS SNSSCITWEGTNGEFKMTDPDEVARRWGERKSKPNMNYDKLSRALRYYYDKNI MTKVH GKRYAYKFDH GIAQALQPH PESSLYKYP SDLPYMG SYHAHPQKMNFV APHPALPVTSSSFFAAPNPYWNSPTGGIYPNTRLPTSHMPSHLGTYY
Purity	Greater than 85% as determined by SDS-PAGE.
Source	E.coli
MW	57.8 kDa
Uniprot ID	P11308
Expiration Date	6 months from date of receipt.