

Product Datasheet

Fungi ACC1 protein (orb604488)

Description	Recombinant Candida albicans Acetyl-CoA carboxylase(ACC1),partial
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged
Note	For research use only
Application notes	Partial
Protein Sequence	GSAMRSLGDKISSTIVAQHAQVPCIPWSGTGVDEVKIDPQTNLVSVADDIYAKGC CTSPEDGLEKAKKIGFPVMIKASEGGGGKGIRKVDDEKNFITLYNQAAANEIPGSPI FIMKLAGDARHLEVQLLADQYGTNISLFGSDCSVQRRHQKIIIEAPVTIARKETFH EMENAAVRLGKLVGYVSAGTVEYLYSHAEDKFYFLELNPRQLQVEHPTTEMVTGV NLPAALQIAMGIPMHRIRDITLYGADPHTTTTDIDFEFKSETSLVSQRRPTPKGH CTACRITSEDPGEGFKPSGGSLHELNFRSSSNVWGYFSVGNQSSIHSFSDSQFG HIFAFGENRQASRKHMVVALKELSIRGDFRTTVEYLIKLETPDFEDNTITTGWLD ELITKKLTAERPDPPIVAVVCGAVTKVHIQAEKEKEYIQSLEKGQVPHRNLLKTIFP VEFIYEGERYKFTATKSSSEDKYTLFLNGSRCVVGARSLSDGGLLCALDGKSHSVY WKEEASATRLSVDGKTCLEVENPTQLRTPSPGKLVKYLVDSGEHVDAGQPYA EVEVMK
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	67.3 kDa
Uniprot ID	A0A1D8PRR7
Expiration Date	6 months from date of receipt.