

Product Datasheet

PLC gamma 1 (phospho-Ser1248) antibody (orb6045)

Description	PLC gamma 1 (phospho-Ser1248) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PLC gamma 1 around the phosphorylation site of Ser1248 (EG(p-S)FE)
Target	Phospho-PLCG1 (Ser1248)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	148 kDa
Uniprot ID	P19174
Expiration Date	12 months from date of receipt.