

Product Datasheet

Bacteria Trichosanthos kirilowii Ribosome-inactivating protein alpha-trichosanthin protein (orb604716)

Description	Recombinant Enterobacteria phage T7 DNA primase/helicase(4),partial
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Partial
Protein Sequence	KKIVVTEGEIDMLTVMELQDCKYPVVS LGHGASAAKKTCAANYEYFDQFEQIILM FDMDEAGRKA VEEAAQVLPAGKVRVAVLPCKDANECHLN GH DREIMEQVWNA GPWIPDGVVSALSLRERIREHLSSEESVGLLFSGCTGINDKTLGARGGEVIMVTS GSGMGKSTFVRQQALQWGTAMGKKVGLAMLEESVEETAEDLIGLHNRVRLRQ SDSLKREIENGKFDQWDFELFGNDTFHLYDSFAEAETDRLLAKLAYMRSGLGC DVIILDHISIVVSASGESDERKMIDNLM TKLGFAKSTGVVLVVICH LKNPDKGKA HEEGRPV SITDLRGSGALRQLSDTI ALERNQQGDMPNLVLVRILKCRFTGDTGIA GYMEYNKETGWLEPSSYS
Purity	Greater than 85% as determined by SDS-PAGE.
Source	E.coli
MW	51.3 kDa
Uniprot ID	P03692
Expiration Date	6 months from date of receipt.