

Product Datasheet

Mouse Adamts13 protein (orb604991)

Description	Recombinant Mouse A disintegrin and metalloproteinase with thrombospondin motifs 13(Adamts13),partia...
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Partial
Protein Sequence	WTPLVGLCSISCGRGLKELYFLCMDSVLKMPVQEELCGLASKPPSRWEVCRARP CPARWETQVLAPCPVTCGGGRVPLSVRCVQLDRGHPISVPHSKCSPVPKPGSFE DCSPEPCPARWKVLSLGPSCASCGLTATQMVACMQLDQGHNEVNETFCKA LVRPQASVPCLIADCAFRWHISAWTECSVSCGDGIQRRHDTCLGPQAQVPVPAN FCQHLPKPMTVRGCWAGPCA
Purity	Greater than 85% as determined by SDS-PAGE.
Source	E.coli
MW	32.7 kDa
Uniprot ID	Q769J6
Expiration Date	6 months from date of receipt.