

Product Datasheet

Mouse Sept9 protein (orb605016)

Description	Recombinant Mouse Septin-9(Sept9)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Full Length
Protein Sequence	MKKSYSYSGVTRTSSGRLRRLADPTGPALKRSFEVEEIEPPNSTPPRRVQTPLLRA TV ASSSQKFQDLGVKNSEPAARLVDSLQSPKPSLRRVELAGAKAPEPMSRRTEIS IDISSKQVESTASAAGPSRFGKRAEVLGHKTPEPVPRRTEITIVKPQESVLRVET PASKIPEGSAPVATDAAPKRVEIQVPKPAEAPNCPLPSQTLSENSEAPMSQLQSRLE PRPSVAEVPYRNQEDSEVTPSCVGDMA DNPRDAMLKQAPASRNEKAPMEFGYV GIDSILEQMRRKAMKQGFENIMVVGQSGLGKSTLINTLFKSKISRKSVQPTSEER IPKTIEIKSITHDIEEKGV RMKLTVIDTPGFGDHINNENCWQPIMKFINDQYEKYLQ EEVNINRKKRIPDTRVHCCLYFIPATGHSRLPLDIEFMKRLSKVVNIVPIAKADTLT LEERVYFKQRITADLLSNGIDVYPQKEFDEDAEDRLVNEKFREMIPFAVVGSDHE YQVNGKRILGRKTKWGTIEVENTHCEFAYLRDLLIRTHMQNIKDITSNIHFEAYR VKRLNEGNSAMANGIEKEPEAQEM
Purity	Greater than 85% as determined by SDS-PAGE.
Source	E.coli
MW	69.6 kDa
Uniprot ID	Q80UG5
Expiration Date	6 months from date of receipt.