

Product Datasheet

Mouse Serpina3n protein (orb605049)

Description	Recombinant Mouse Serine protease inhibitor A3N(Serpina3n)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	Full Length of Mature Protein
Protein Sequence	SFPDGTGLGMDAAVQEDHDNGTQLDSLTLASINTDFAFSLYKELVLKNPDKNIVFS PLSISAALAVMSLGAKGNTLEEILEGLKFNLTTETSEADIHQGFHLLQRLNQPKDQ VQISTGSALFIEKRQQILTEFQEKAKTLYQAEFTADFQQPRQAKKLINDYVRKQT QGMKELVSDLDKRTLMLVNVYIFKAKWKVPFDPLDTFKSEFYAGKRRPVIVPM MSMEDLTTPYFRDEELSCTVVELKYTGNASALFILPDQGRMQQVEASLQPETLRK WKNSLKPRMIDELHLPKFISTDYSLEDVLSKLGIREVFSTQADLSAITGTDLRVS QVVHKAVLDVAETGTEAAAATGVKFPMSAKLYPLTVYFNRPFLLIMIFDTETEIAP FIAKIANPK
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	60.8 kDa
Uniprot ID	Q91WP6
Expiration Date	6 months from date of receipt.