

Product Datasheet

Human PELI3 protein (orb605072)

Description	Recombinant Human E3 ubiquitin-protein ligase pellino homolog 3 (PELI3)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Full Length
Protein Sequence	MVLEGNPEVGSPRTSDLQHRGNKGSCVLSSPGEDAQPGEPIKYGELIVLGCCE EGGEETEAQRGEVTGPRAHSCYNGCLASGDKGRRRSRLALSRRSHANGVKPDV MHHISTPLVSKALSNRGQHSISYTLRSHSVIVEYTHSDTDMFQIGRSTENMIDF VVTDTSPGGGAAEGPSAQSTISRYACRILCDRPPYTARIYAAGFDASSNIFLGER AAKWRTPDGLMDGLTTNGVLVMPAGGFSEDSAPGVWREISVCGNVYTLRDS RSAQQRGKLVENESNVLDQDGLDLGATLLWRTPAGLLRAPTLKQLEAQRQEA NAARPQCPVGLSTLAFPSARGRTAPDKQQPWVYVRCGHVHGYHGWGCRER GPQERECPLCRLVGPYVPLWLGQEAGLCLDPGPPSHAFAPCGHVCSEKTARYW AQTPLPHGTHAFHAACPFCAWLTGEHGCVRILIFQGPLD
Purity	Greater than 85% as determined by SDS-PAGE.
Source	E.coli
MW	54.8 kDa
Uniprot ID	Q8N2H9
Expiration Date	6 months from date of receipt.