

Product Datasheet

Human PELI1 protein (orb605083)

Description	Recombinant Human E3 ubiquitin-protein ligase pellino homolog 1 (PELI1)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Full Length
Protein Sequence	MFSPDQENHPSKAPVKYGEIVLGYNGSLPNGDRGRRKSRFALFKRPKANGVKP STVHIACTPQAAKAISNKDQHSISYTLSTRAQTVVVEYTHDSNTDMFQIGRSTESPI DFVVTDTVPGSQNSDTQSVQSTISRFAICRIICERNPPFTARIYAAGFDSSKNIFL GEKAAKWKTSQMDGLTTNGVLVMHPRNGFTEDSKPGIWREISVCGNVFSLR ETRSAQQRGKMVEIETNQLQDGLDLCGATLLWRTAEGLSHTPTVKHLEALRQ EINAARPQCPVGFNTLAFPSMKRKDVVDEKQPWVYLNCGHVHGYHNWGNKEE RDGKDRECPMCRSVGPYVPLWLGCEAGFYVDAGPPTHAFSPCGHVCSEKTTAY WSQIPLPHGTHTFHAACPFAHQQLAGEQGYIRLIFQGPLD
Purity	Greater than 85% as determined by SDS-PAGE.
Source	E.coli
MW	50.4 kDa
Uniprot ID	Q96FA3
Expiration Date	6 months from date of receipt.