

Product Datasheet

Rat Cfi protein (orb605203)

Description	Recombinant Rat Complement factor I(Cfi)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Full Length of Mature Protein
Protein Sequence	KNTPASGQPQEDLVEQKCLLNKYNTHHSCDKVFCQPWQKCIEGTCACKLPYQCPK AGTPVCATNGRGYPTYCHLKSFECLHPEIKFSNNGTCTAEEKFNVS LIYGSTDTE GIVQVKLVDQDEKMFICKNSWSTVEANVACFDLGFPLGVRDIQGRFNIPVNHKIN STECLHVRCQGVETSLAECTFTKKSSKAPHGLAGVVCYTQDADFPTSQSFCVN GKRIPQEACDGVND CGDQSDCLCKGCRGQAFLCKSGVCIPNQRKCNGEVD CITGEDES GCEEDKKNKIHKGLARSDQGGETEIEETEEMLT PDM DTERKRIKSL LPKLSCGVKRNTHIRRKRVVGGKPAEMGDYPWQVAIKDGD RITCGGIYIGGCWIL TAAHCVRPSRYRNYQVWTSLLDWLKPNSQLAVQGVS R VVHEKYN GATYQNDI ALVEMKKHPGKKECELINSVPACVPWSPYLFQPNDRCIISGWGREKDNQKVYSL RWGEVDLIGNCSRFPGRYYEKEMQCAGTSDGSIDACKGDSGGPLVCKDVNNV TYVWGVSWGENCGKPEFPGVYTRVASYFDWISYYVGRPLVSQYNV
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Mammalian cell
MW	70.4 kDa
Uniprot ID	Q9WUW3
Expiration Date	6 months from date of receipt.