

Product Datasheet

Human CSF1R protein (orb605208)

Description	Recombinant Human Macrophage colony-stimulating factor 1 receptor(CSF1R),partial
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal hFc-tagged
Note	For research use only
Application notes	Partial
Protein Sequence	IPVIEPSVPELVVKPGATVTLRCVNGSVIEWDGPSPHWTLYSDGSSSILSTNNA TFQNTGTYRCTEPGDPLGGSAIHLVYKDPARPWNVLAQEVVVFEDQDALLPCL LTDPVLEAGVSLVRVRGRPLMRHTNYSFSPWHGFTIHRAKFIQSQDYQCSALMG GRKVMISISIRLKVKVQVIPGPPALTLVPAELVRIRGEAAQIVCSASSVDVNFDFVLQ HNNTKLAIPQQSDFHNNRYQKVLTNLNDQVDFQHAGNYSVASNVQGKHSTSM FFRVVESAYLNLSSQNLIQEVTVGEGNLKVMVEAYPGLOGFNWYTLGPFSDH QPEPKLANATTKD TYRHTFTLSLPR LKPSEAGRYSFLARNPGGW RALTFELTLRY PPEVSVIWT FINGSGTLLCAASGYQP NVTWLQCSGHTDRCDEA QVLQVWDDP YPEVLSQEPFHKVTVQSLTVETLEHNQTYECRAHNSVSGSGSWAFIPISAGAHTH PPDE
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Mammalian cell
MW	82.0 kDa
Uniprot ID	P07333
Expiration Date	6 months from date of receipt.