

Product Datasheet

Rat Pkm protein (orb605241)

Description	Recombinant Rat Pyruvate kinase PKM(Pkm)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Full Length
Protein Sequence	MPKPDSEAGTAFIQTLQLHAAMADTFLEHMCRLDIDSAPITARNTGIICTIGPASR SVEMLKEMIKSGMNVARLNFSGHTHEYHAETIKNVRAATESFASDPILYRPVAVA LDTKGPEIRTGLIKSGTAEVELKKGATLKITLDNAYMEKCDENILWLDYKNICKV VEVGSKIYVDDGLISLQVKEKGADYLVTEVENGGSLGSKKGVNLPGAADVLPV SEKDIQDLKFGVEQDVMVFASFIRKAADVHEVRKVLGEKGKNIKISKIENHEGV RRFDEILEASDGIMVARGDLGIEIPAENVFLAQMMIGRCNRAGKPVICATQMLES MIKKPRPTRAEGSDVANAVLDGADCIMLSGETAKGDYPLEAVRMQHIAAREA AVFHRLLEELARASSQSTDPLEAMAMGSVEASYKCLAAALIVLTESGRSAHQVA RYRPRAPIIAVTRNPQTARQAHLRYGIFPVLCCKDAVLDAWAEDVDLRVNLAMNV GKARGFFKKGDVVIVLTGWRPGSGFTNTMRVVPV
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Mammalian cell
MW	61.8 kDa
Uniprot ID	P11980
Expiration Date	6 months from date of receipt.