

Product Datasheet

Rabbit TFPI protein (orb605252)

Description	Recombinant Rabbit Tissue factor pathway inhibitor(TFPI)
Endotoxins	Less than 1.0 EU/ug as determined by LAL method.
Preservatives	Lyophilized from a 0.2 µm filtered 20 mM Tris-HCl, 0.5 M NaCl, 6% Trehalose, pH 8.0
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged
Note	For research use only
Application notes	Full Length of Mature Protein
Protein Sequence	AAEEDDEFTNITDIKPLQKPTHSFCAMKVDDGPCRAYIKRFFFNILTHQCEEFIYG GCEGNENRFESLEECKEKCARDYPKMTTKLTFQKGKPDFCFLEEDPGICRGYITR YFYNNQSKQCERFKYGGCLGNLNNFESLEECKNTCENPTSDFQVDDHRTLNT VNNTLINQPTKAPRRWAFHGPSWCLPPADRGLCQANEIRFFYNAILGKCRPFKYS GCGGNENNFTSKKACITACKKGFIPKSIKGLIKTRKKKKQPVKITYVETVVKKT
Purity	Greater than 95% as determined by SDS-PAGE.
Source	Mammalian cell
MW	35.4 kDa
Uniprot ID	P19761
Expiration Date	6 months from date of receipt.