

Product Datasheet

Mouse Pkdcc protein (orb605275)

Description	Recombinant Mouse Extracellular tyrosine-protein kinase PKDCC(Pkdcc)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Full Length of Mature Protein
Protein Sequence	PRPGQSPGSSAAPGPGRRRGRELARQIRERYEEVQRYSRGGPGGAGRPERR RLMDLAPGGPGLQRP RPVRSPPDGAPGWPPAPGPGSPGPGPRLGCAALRNV SGAQYVSGSYTKAVYRVRLPGGAAVALKAVDFSGHDLGSCVREFGARRGCYRL AAHKLLKEMVLLERLRHPNVLQLYGYCYQDSEGIPDTLTITELGAPVEMIQLLQT SWEDRFRIKLSLGRLLHHLAHSPLGVSVTLLDFRPRQFVLVNGELKVTDLDDARVE ETPCTSSADCTLEFPARNFSLPCSAQGWCEGMNEKRNLNAYRFFFTYLLPHSAP PSLRPLLDIVNATGELAWGVDETLAQLETALHLFRSGQYLQNSTSSRAEYQRIP DSAITQEDYRCWPSYHHGGCLLSVFNLAEAIDVCESHAQCRAFFVTNQTWTG RKLVFVKTGWNQVVPDAGKTTYVKAPG
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Mammalian cell
MW	55.5 kDa
Uniprot ID	Q5RJI4
Expiration Date	6 months from date of receipt.