

Product Datasheet

Virus polA protein (orb605289)

Description	Recombinant Thermus aquaticus DNA polymerase I, thermostable(polA),partial
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Application notes	Partial
Protein Sequence	ESPKALEEAPWPPPEGAFVGFVLSRKEPMWADLLALAAARGGRVHRAPEPYKAL RDLKEARGLLAKDLSVLALREGLGLPPGDDPMLLAYLLDPSNTTPEGVARRYGG EWTEEAGERAALSERLFANLWGRLEGEERLLWLYREVERPLSAVLAHMEATGV RLDVAYLRALSLEVAEEIARLEAEVFRLAGHPFNLNSRDQLERVLFDLGLPAIGK TEKTGKRSTSAAVLEALREAHPIVEKILQYRELTKLKSTYIDPLPDLIHPRTGRLHTR FNQTATATGRLSSSDPNLQNIPVRTPLGQRIRRAFIAEEGWLLVALDYSQIELRVL AHLSGDENLIRVFQEGRDIHTETASWMFGVPREAVDPLMRRAAKTINFGVLYGM SAHRLSQELAIPIYEEAQAFIERYFQSFPKVRWIEKTLLEGRRRGYVETLFGRRRY VPDLEARVKSVEREAERMAFNMPVQGTAAADLMKLAMVKLFPRLEEMGARMLLQ VHDELVLAPKERAEEAVARLAKEVMGVYPLAVPLEVEVGIGEDWLSAKE
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	88.3 kDa
Uniprot ID	P19821
Expiration Date	6 months from date of receipt.