

Product Datasheet

Rat Nefl protein (orb605354)

Description	Recombinant Rat Neurofilament light polypeptide(Nefl)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Full Length of Mature Protein
Protein Sequence	SSFSYEPYFSTSYKRRYVETPRVHISSVRSGYSTARSAYSYSAPVSSSLSVRRSY SSSSGSLMPSENLDLSQVAAISNDLKSIRTQEKAQLQDLNDRFASFIERVHELE QQNKVLEAELLVLRQKHSEPSRFRALYEQEIRDLRLAAEDATNEKQALQGERGL EETLRNLQARYEEEVLSREDAEGRLEARKGADEAALARAEELEKRIDSLMDEIAF LKKVHEEEIAELQAQIQYAQISVEMDVSSKPDLSAALKDIRAQYEKLAAKNMQNA EEWFKSRFTVLTESAAKNTDAVRAAKDEVSESRLLKAKTLEIEACRGMNEALEK QLQELEDKQNADISAMQDTINKLENELRSTKSEMARYLKEYQDLLNVKMALDIEI AAYRKLLEGEETRLSFTSVGSITSQSSQVFGRSAYSGLSQSSSYLMSARAFPA YYTSHVQEEQSEVEETIEATKAEEAKDEPPSEGEAEKEKEKEKEKEKEKEKEKEKE AAKDESEDAKEEEGGEGEEEDTKESEEEKEKEESAGEEQAAKKKD
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Yeast
MW	63.2 kDa
Uniprot ID	P19527
Expiration Date	6 months from date of receipt.