

Product Datasheet

Human MADCAM1 protein (orb605485)

Description	Recombinant Human Mucosal addressin cell adhesion molecule 1(MADCAM1),partial
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal Myc-tagged and C-terminal 6xHis-tagged
Note	For research use only
Application notes	Partial
Protein Sequence	QSLQVKPLQVEPPEPVVAVALGASRQLTCRLACADRGASVQWRGLDTSLGAVQ SDTGRSVLTVRNASLSAAGTRVCVGSCGGRTFQHTVQLLVYAFPDQLTVSPAAL VPGDPEVACTAHKVTPVDPNALSFSLLVGGQELEGAQALGPEVQEEEEEPQGDE DVLFRVTERWRLPPLGTPVPPALYCQATMRLPGLELSHRQAIPVLHSPTSPEPPD TTSPEPDTTSPESPDTTSQEPPDTSPEPPDKTSPEPAPQQGSTHTPRSPGSTR TRRPEISQAGPTQGEVIPTGSSKPAGDQ
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	34.4 kDa
Uniprot ID	Q13477
Expiration Date	6 months from date of receipt.