

Product Datasheet

Plant UGT89C1 protein (orb605533)

Description	Recombinant Arabidopsis thaliana UDP-glycosyltransferase 89C1(UGT89C1)
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged
Note	For research use only
Application notes	Full length
Protein Sequence	MTTTTTKKPHVLVIPFPQSGHMPHLDLTHQILLRGATVTVLVTPKNSSYLDALRS LHSPHFKTLILPFPSPHCIPSGVESLQQLPLEAIVHMFDAISRHLHDPLVDFLSRQP PSDLPDAILGSSFLSPWINKVADAFSISISFLPINAHISVMMWAQEDRSFFNDLET ATTESYGLVINSFYDLEPEFVETVKTRFLNHHRIWTVGPLLPFKAGVDRGGQSSIP PAKVSALWDSCPEDNSVYVVGFGSQIRLTAEQTAALAAALEKSSVRFIWAVERDA AKKVNSSDNSVEEDVIPAGFEERVKEKGLVIRGWAPQTMILEHRAVGSYLTHLG WGSVLEGMVGGVMMLLAWPMQADHFFNTTLIVDKLRAAVRVGENRDSVPDSDK LARILAESAREDLPERVTLMKLREKAMEAIKEGGSSYKNLDELVAEMCL
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Yeast
MW	50.6 kDa
Uniprot ID	Q9LNE6
Expiration Date	6 months from date of receipt.