

Product Datasheet

FOXO4 (phospho-Ser197) antibody (orb6056)

Description	FOXO4 (phospho-Ser197) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human FOXO4 around the phosphorylation site of Ser197 (AA(p-S)MD)
Target	phospho-FOXO4 (Ser197)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54 kDa
Uniprot ID	P98177
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.