

Product Datasheet

p63 Antibody / Tumor protein 63 (orb606286)

Description

p63 is a homolog of the tumor suppressor

epithelial laye...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

A portion of amino acids 3-106 were used as the immunogen for this p63 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the p63 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The optimal dilution of the p63 antibody for each application should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2b, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

Q9H3D4

Hazard

Information

This p63 antibody is available for research use only.

Dilution Range

ELISA: order BSA/sodium azide-free format for coating, Western blot: 1-2ug/ml

Expiration Date

12 months from date of receipt.