

## Product Datasheet

**TRAcP Antibody / TRAP / ACP5 (orb606365)**

## Description

It recognizes a protein of 35kDa, which is

phosphatase (TRAcP)...

## Species/Host

Rabbit

## Reactivity

Human, Mouse, Rat

## Conjugation

Unconjugated

## Tested

IHC-P

## Applications

## Immunogen

Recombinant full-length human protein was used as the immunogen for the recombinant TRAcP antibody.

## Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

## Storage

Store the recombinant TRAcP antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

## Note

For research use only

## Application notes

Optimal dilution of the recombinant TRAcP antibody should be determined by the researcher.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

## Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

## Isotype

Rabbit IgG, kappa

## Clonality

Recombinant

## Purity

Protein A affinity chromatography

## Uniprot ID

**P13686**

## Hazard

## Information

This recombinant TRAcP antibody is available for research use only.

## Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

## Expiration Date

12 months from date of receipt.