

Product Datasheet

Keratin 6 Antibody / Cytokeratin 6 (orb606593)

Description

This mAb recognizes a protein of 56kDa,

humans, multiple isofo...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FACS, IHC-P

Applications
Immunogen

Recombinant human KRT6A protein was used as the immunogen for the Cytokeratin 6 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Cytokeratin 6 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the Cytokeratin 6 antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2a, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P02538

Hazard
Information

This Cytokeratin 6 antibody is available for research use only.

Dilution Range

Flow cytometry: 1-2ug/million cells, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.