

## Product Datasheet

### p63 Antibody / Tumor protein 63 (orb606933)

## Description

p63 is a homolog of the tumor suppressor

epithelial laye...

## Species/Host

Mouse

## Reactivity

Human, Mouse

## Conjugation

Unconjugated

## Tested

IHC-P

## Applications

## Immunogen

Full length human recombinant protein was used as the immunogen for the p63 antibody.

## Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

## Storage

Store the p63 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

## Note

For research use only

## Application notes

Optimal dilution of the p63 antibody should be determined by the researcher.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

## Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

## Isotype

Mouse IgG2a, kappa

## Clonality

Monoclonal

## Purity

Protein G affinity chromatography

## Uniprot ID

**Q9H3D4**

## Hazard

## Information

This p63 antibody is available for research use only.

## Dilution Range

Immunohistochemistry (FFPE): 0.25-0.5ug/ml for 30 min at RT

## Expiration Date

12 months from date of receipt.