

Product Datasheet

ENAH Antibody / Mena (orb606938)

Description

The Wiskott-Aldrich syndrome (WAS) is

defects in cell-med...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

A portion of amino acids 485-589 from the human protein was used as the immunogen for the ENAH antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the ENAH antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The stated application concentrations are suggested starting points. Titration of the ENAH antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG2c, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

Q8N8S7

Hazard

Information

This ENAH antibody is available for research use only.

Dilution Range

ELISA (order BSA/sodium azide-free format for coating), Western blot: 1-2ug/ml

Expiration Date

12 months from date of receipt.