

Product Datasheet

Bcl2L2 Antibody / Bcl-W (orb606954)

Description

Bcl2L2, also called BCLW, promotes cell

apoptosis. Mediates s...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

Full-length recombinant human protein was used as the immunogen for this Bcl2L2 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the Bcl2L2 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the Bcl2L2 antibody to be titrated up or down for optimal performance.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

[Q92843](#)

Hazard

Information

This Bcl2L2 antibody is available for research use only.

Dilution Range

Western blot: 1-2ug/ml

Expiration Date

12 months from date of receipt.