

Product Datasheet

**PD-L2 Antibody / Programmed death ligand 2 /
CD273 (orb607156)**

Description

Recognizes a protein identified as PD-

of CD28 by B7-1 (CD80) or B7...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FACS, IF, WB

Applications
Immunogen

A recombinant fragment of human PD-L2 protein within amino acids 27-220 was used as the immunogen for the PD-L2 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the PD-L2 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the PD-L2 antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

Q9BQ51

Hazard
Information

This PD-L2 antibody is available for research use only.

Dilution Range

ELISA (order BSA/sodium azide-free format for coating), Flow cytometry: 1-2ug/million cells, Immunofluorescence: 1-2ug/ml, Western blot: 1-2ug/ml

Expiration Date

12 months from date of receipt.