

Product Datasheet

CD66 antibody (CF488) (orb607216)

Description

Mouse monoclonal antibody against

conjugated to CF488

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF

Applications
Immunogen

Recombinant full-length human CEA protein

Target

CEACAM1 & CEACAM5

Preservatives

Antibody Purified from Bioreactor Concentrate by Protein A/G and conjugated to various reporter molecules. Prepared in 10mM PBS with 0.05% rAlbumin and 0.05% azide. Contact us if you require this Ab in a different format.

Concentration

Purified Ab conjugated to CF488

Storage

Antibody with azide - store at 4 to 8°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Mouse / IgG2a, kappa

Clonality

Monoclonal

MW

80-200kDa

Uniprot ID

P06731

Entrez

634, 1048

Dilution Range

Flow Cytometry (5ul per test per one million cells or 5ul per 100ul of whole blood);Immunofluorescence (1:50-1:100);Optimal dilution for a specific application should be determined.

Expiration Date

12 months from date of receipt.