

Product Datasheet

Ku (phospho-70/80) antibody (CF488) (orb607366)

Description

Mouse monoclonal antibody against Ku

	CF488
Species/Host	Mouse
Reactivity	Human, Primate, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF
Immunogen	Nuclear extract of human HL-60 cells
Target	XRCC6 (p70) & XRCC5 (p80)
Preservatives	Antibody Purified from Bioreactor Concentrate by Protein A/G and conjugated to various reporter molecules. Prepared in 10mM PBS with 0.05% rAlbumin and 0.05% azide. Contact us if you require this Ab in a different format.
Concentration	Purified Ab conjugated to CF488
Storage	Antibody with azide - store at 4 to 8°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	70kDa & 80kDa
Uniprot ID	P13010, P12956
Entrez	7520, 2547
Dilution Range	Flow Cytometry (5ul per test per one million cells or 5ul per 100ul of whole blood);Immunofluorescence (1:50-1:100);Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.