

Product Datasheet

GLG1 antibody (CF700) (orb607386)

Description

Mouse monoclonal antibody against GLG1

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, IF

Immunogen

Golgi fraction from human liver cells

Target

GLG1

Preservatives

Antibody Purified from Bioreactor Concentrate by Protein A/G and conjugated to various reporter molecules. Prepared in 10mM PBS with 0.05% rAlbumin and 0.05% azide. Contact us if you require this Ab in a different format.

Concentration

Purified Ab conjugated to CF700

Storage

Antibody with azide - store at 4 to 8°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

134kDa

Uniprot ID

[Q92896](#)

Entrez

[2734](#)

Dilution Range

Flow Cytometry (5ul per test per one million cells or 5ul per 100ul of whole blood); Immunofluorescence (1:50-1:100)Optimal dilution for a specific application should be determined.

Expiration Date

12 months from date of receipt.