

Product Datasheet

KRT18 antibody (orb607502)

Description

Mouse monoclonal antibody to Cytokeratin 18

Species/Host	Mouse
Reactivity	Bovine, Canine, Equine, Feline, Ferret, Goat, Hamster, Human, Monkey, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Cytoskeletal preparation extracted from the human vulvar squamous carcinoma cell line A431
Target	KRT18
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	45kDa
Uniprot ID	P05783
Entrez	3875
Dilution Range	Flow Cytometry (1-2ug/million cells); Immunofluorescence (1-2ug/ml); Western Blot (1-2ug/ml); ,Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 min at RT),(Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes); Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709. United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558