

Product Datasheet

TACSTD1 antibody (PE) (orb607511)

Description

Mouse monoclonal antibody against Ep-CAM /

Marker) conjugat...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, IF

Immunogen

Small cell lung carcinoma cells

Target

TACSTD1

Preservatives

Antibody Purified from Bioreactor Concentrate by Protein A/G and conjugated to various reporter molecules. Prepared in 10mM PBS with 0.05% rAlbumin and 0.05% azide. Contact us if you require this Ab in a different format.

Concentration

Purified Ab conjugated to PE

Storage

Antibody with azide - store at 4 to 8°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

40-43kDa

Uniprot ID

P16422

Entrez

4072

Dilution Range

Flow Cytometry: 5ul per test per one million cells (or 5ul per 100ul of whole blood);Immunofluorescence (1:50-1:100)Optimal dilution for a specific application should be determined.

Expiration Date

12 months from date of receipt.