

Product Datasheet

TP53 antibody (orb607761)

Description

Mouse monoclonal antibody to p53 Tumor

Species/Host	Mouse
Reactivity	Bovine, Canine, Hamster, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Gel-Purified p53-beta-galactosidase fusion protein containing murine p53 from aa 14-389
Target	TP53
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	53kDa
Uniprot ID	P04637
Entrez	7157
Dilution Range	ELISA (For coating, order Ab without BSA)Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.