

Product Datasheet

VEGFA antibody (orb607789)

Description

Mouse monoclonal antibody to VEGF (Vascular

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Human VEGF189 recombinant protein
Target	VEGFA
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	19-22kDa (reducing) and 38kDa-44kDa (non-reducing)
Uniprot ID	P15692
Entrez	7422
Dilution Range	ELISA (For coating, order Ab without BSA); Immunofluorescence (1-2ug/ml); Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.