

Product Datasheet

PDCD1LG2 antibody (orb607813)

Description

Mouse monoclonal antibody to PD-L2 /

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Immunogen	Recombinant fragment (around aa 27-220) of human PD-L2 protein (exact sequence is proprietary)
Target	PDCD1LG2
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	31kDa
Uniprot ID	Q9Q51
Entrez	80380
Dilution Range	ELISA (For coating, order antibody without BSA);,Flow Cytometry (1-2ug/million cells);,Immunofluorescence (1-2ug/ml);,Western Blot (1-2ug/ml),Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.