

## Product Datasheet

### TNFRSF10B antibody (orb607824)

**Description**

Mouse monoclonal antibody to CD262 / DR5

**Species/Host**

Mouse

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

IHC-P

**Applications**
**Immunogen**

Recombinant human DR5 protein fragment (around aa266-393) (exact sequence is proprietary)

**Target**

TNFRSF10B

**Preservatives**

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

**Concentration**

Purified Ab with BSA and Azide at 200ug/ml

**Storage**

Antibody with azide - store at 2 to 8°C.  
Antibody without azide - store at -20 to -80°C.  
Antibody is stable for 24 months. Non-hazardous. No MSDS required.

**Note**

For research use only

**Isotype**

Mouse / IgG1, kappa

**Clonality**

Monoclonal

**MW**

48kDa

**Uniprot ID**

521456

**Entrez**

**8795**

**Dilution Range**

Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 min at RT),(Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes); Optimal dilution for a specific application should be determined.

**Expiration Date**

12 months from date of receipt.