

Product Datasheet

MS4A1 antibody (CF488) (orb607866)

Description

Mouse monoclonal antibody against CD20 /

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Lymphoblastoid cell line Daudi
Target	MS4A1
Preservatives	Antibody Purified from Bioreactor Concentrate by Protein A/G and conjugated to various reporter molecules. Prepared in 10mM PBS with 0.05% rAlbumin and 0.05% azide. Contact us if you require this Ab in a different format.
Concentration	Purified Ab Conjugated to CF488
Storage	Antibody with azide - store at 4 to 8°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG2a, kappa
Clonality	Monoclonal
MW	33-37kDa
Uniprot ID	P11836
Entrez	931
Dilution Range	Flow Cytometry: 5ul per test per one million cells (or 5ul per 100ul of whole blood)Immunofluorescence (1:50-1:100)Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.